



General Wiring Notes

EASYLAB Controller TCU3



Last update
16.02.2024

Notes

The options contained in these wiring instructions can be activated and parameterized with the EasyConnect software. It is not possible to set the controller parameters individually at the factory.

Applicable documents

Leaflets and assembly instructions of

- EASYLAB Controller TCU3
- EASYLAB Expansion modules EM-TRF or EM-TRF-USV, EM-LON, EM-AUTOZERO, EM-LIGHT, EM-BAC-MOD-01, EM-IP
- Sensors/Transducers for fume cupboard control VS-TRD or DS-TRD-02
- EASYLAB Control panels BE-SEG-02, BE-LCD, CP-TOUCH

General safety instructions

Recognised codes of engineering practice, especially the safety requirements and accident prevention regulations must be observed when performing all assembly, wiring and commissioning procedures.

Only specially trained and qualified personnel are permitted to perform the assembly, wiring and commissioning. The electrical installations must be implemented in compliance with the current valid EN / VDE / DIN and local electrical regulations.

Danger:

Electric shock by touching live parts.

Electrical devices may cause electrical hazards during operation

Before installation/assembly works begin:

**Disconnect all poles or switch off supply voltage of EASYLAB controller.
Protect power supply against unauthorised or unintended restart.**

All installation operations may only start after these safety operations.





Contents




Overview terminal assignment	3
Power supply	4
EASYLAB Communication line (CL) - sockets	5
EASYLAB Communication line (CL) - screw terminals	6
EASYLAB Communication system structure and network termination	7
EASYLAB system design and Room Management Function (RMF)	8
EASYLAB Control panels	9
Service socket for configuration and diagnostics	10
Digital inputs (DI) / Digital outputs (DO)	11
Analogue inputs (AI) / Analogue outputs (AO)	12
Analogue outputs (AO) / Terminal 3 (connection TROX HPD) / DS-TRD	13
Input / Output assignment of Fume Cupboard Controller (FH)	14
Input / Output assignment of supply air / extract air Controller (SC/EC)	15
Input / Output assignment of Room Controller (Supply/Extract air)	16
Input / Output assignment of TROX Adaption Module (TAM)	17
Wiring schematic Fume Cupboard Controller with face velocity transducer (FH-VS)	18
Wiring schematic Fume Cupboard Controller with sash distance sensor (FH-DS, FH-DV)	19
Wiring schematic Fume Cupboard Controller with face velocity transducer & sash distance sensor (FH-VD)	20
Wiring schematic Fume Cupboard Controller with 3-Point control (FH-3P)	21
Wiring schematic Fume Cupboard Controller with 2-Point control (FH-2P)	22
Wiring schematic Fume Cupboard Controller with fixed value control (FH-F)	23
Fume Cupboard controller special functions Part 1 (Sash monitoring, Fume scrubber, Supportive flow technology, sash close device)	24
Fume Cupboard controller special functions Part 2 (Mode dependent switching, Fume cupboard lightning, Smoke extraction)	25
Fume Cupboard controller special functions Part 3 (Volume flow integration, Temperature alarm, Motion sensor)	26
Connection diagram Supply air-/ extract air controller 0-10 V / 2-10 V (SC-E0; SC-E2; EC-E0; EC-E2)	27
Connection diagram Supply air-/ extract air controller with 3-point-regulation (SC-3P; EC-3P)	28
Connection diagram Supply air-/ extract air controller with 2-point-regulation (SC-2P; EC-2P)	29
Connection diagram Supply air-/ extract air controller with fixed volume flow (SC-F; EC-F)	30
Room controller special functions Part 1 (Sun blind control, Lightning control, Supply air diffusor control, Volume flow integration)	31
Room controller special functions Part 2 (Operation mode dependent switching)	32
EASYLAB Pressure control	33
Pressure control connections; tubing	34
LonWorks [®] Interface (Expansion module EM-LON)	35
BACnet MS/TP, Modbus RTU Interface (Expansion module EM-BAC-MOD-01)	36
BACnet IP, Modbus IP Interface (Expansion module EM-IP)	37




Exterior sockets / Indicators:


Left box side



Right box side



Interior sockets / Indicators:





- 1 Red error indicator
- 2 Yellow indicator - Network termination activated
- 3 Green indicator (Reserved)
- 4 Yellow indicator CL-Network data received
- 5 Green indicator Controller operation (Heartbeat)
- 6 Socket X1 DI1 - Sash contact 500 mm
- 7 Socket X2 Control panel 1
- 8 Socket X3 Control panel 2
- 9 Socket X4 A04 Actuator
- 10 Socket X5 AI5 - Face velocity transducer
- 11 Socket X6 Communication line 1 - RJ45
- 12 Socket X7 Communication line 2 - RJ45

- 13 Socket Fume cupboard lights via connector from Wago Winsto Midi Series 770-113 (only with EM-LIGHT expansion module)
- 14 Socket Magnetic valve EM-AUTOZERO
- 15 Socket AI1 for internal volume flow transducer
- 16 Expansion slot for LonWorks®, BACnet, Modbus interface with EM-LON, EM-BAC-MOD-01 or EM-IP
- 17 Optical status indicators of digital inputs
- 18 Activation switch for CL network termination
- 19 Screw terminals for communication line 1
- 20 Screw terminals for communication line 2

- 21 PE connection
- 22 Clamp / Connection pad for cable shield
- 23 Screw terminals Analogue inputs AI1...AI4
- 24 connection digital actuator
- 25 Screw terminals temperature sensor PT1000 T11
- 26 Screw terminals Digital inputs DI2...DI6
- 27 Screw terminals Analogue outputs AO1...AO3
- 28 Screw terminals Power supply 24 V AC / DC
- 29 Screw terminals Digital outputs DO1... DO6
- 30 Green indicator 24 V ok
- 31 Connection pad for Expansion module mains supply
- 32 Optical indicators for switching state of digital outputs




Power supply



24 V Power LED: On = Supply voltage ok
Off = Supply voltage not ok

Supply voltage 24 V AC / DC



Technical data:

- Supply voltage 24 V AC \pm 15% 50-60 Hz
- Supply voltage 24 V DC \pm 15%
- Fuse 2,5 A slow, 250 V; glass tube fuse 5x20 mm
- Power consumption up to 40 VA (Maximum configuration) inrush current actuator: I_{max} 20 A for 5 ms
- Only TCU3 without additional configuration:5 VA
- Maximum power consumption is determined by the individual controller configuration. Typical configurations will lead to following values:
- Fume cupboard controller with control panel up to 35 VA
- Room controller with room control panel CP-TOUCH..... up to 33 VA
- Room controller without room control panel up to 29 VA
- Supply air/extract air controller without room control panel up to 29 VA

Important wiring instructions:

Observe polarity

Polarity of DC and AC supply voltage has to be strictly complied for all controllers during wiring.

Limited wiring in a row

Power supply wiring in a row (via double terminals) is only approved for a maximum of 5 EASYLAB controllers!

Supply voltage at analogue inputs for supply of sensor/transducer


Type of used controller power supply (AC/DC) will influence type of provided power supply at analogue inputs:

- TCU3 Power supply 24 V AC \rightarrow Provided sensor supply at AI1-AI5 is 24 V AC
- TCU3 Power supply 24 V DC \rightarrow Provided sensor supply at AI1-AI5 is 24 V DC

This may require attention for connected sensors

Mains supply voltage 230 V AC

Alternative power supply
(only with expansion module EM-TRF or EM-TRF-USV)



Technical data:

- Supply voltage 230 V AC \pm 20% 50-60 Hz
- Fuse 500 mA slow, 250 V
- Power consumption up to 40 VA (Maximum configuration)

Important wiring instruction:

- Never connect 230 V and 24 V supply voltage simultaneous at one controller.
- Never connect 24 V supply voltage at controller with built in expansion module EM-TRF or EM-TRF-USV

Supply voltage at analogue inputs for supply of sensor/transducer

TCU3 Power supply 230V AC \rightarrow Provided sensor supply at AI1-AI5 is 24 V DC (with EM-TRF, EM-TRF-USV)

This may require attention for connected sensors



Communication line (CL)

For Plug- & Play data exchange between the controllers within one EASYPAB system



Technical data:

- Maximum device no. within one system: 24
- Maximum total length: 300 m
- Network cable type **SF-UTP** (braid + foil shield) according ISO IEC 11801 (2002)
patch cable with RJ45 plug connectors on both sides or cable from a reel; e.g. Cat 5
- Activate termination resistor at the beginning and the end of the EASYPAB communication line CL.
- 5 m connection cable for communication line CL included in the delivery
(only for fume cupboard controller with device function FH-xxx and TAM)

Note:


Use always green patch cables for the communication line CL

Aim:

Simple determination of communication line cable (green) and control panel connecting cable (blue)



CL - Connection sockets and indicators



Termination indicator

Yellow = Termination active Off = Termination inactive

Data reception CL

Yellow / blinking = Data received

Heartbeat

Green slow blinking = Normal controller operation
Green off = Device not operable

Default socket communication line 1

Default socket communication line 2

CL - Termination resistor



Activation switch
CL termination

Termination of communication line CL

For a correct data exchange of the controllers is a termination on both ends of the communication line required.



ON = Termination activated
OFF = Termination deactivated





Communication line (CL)

CL - Screw terminals as an alternative wiring option:

For data exchange via network cable from a reel as an alternative to patch cables


Important wiring note:



For port COMM1 or COMM2 on one EASYPAB controller use either the RJ45 sockets with patch cable **or** the screw terminals with network cable from a reel to built up the communication line CL.

Technical data:

- Maximum device no. within one system: 24
- Maximum total length: 300 m
- Network cable type **SF-UTP** (braid + foil shield) according ISO IEC 11801 (2002) cable on a reel; e.g. Cat 5
- Activate termination resistor at the beginning and the end of the EASYPAB communication line CL.




Screw terminals COMM1
COMM-H / COMM-L / GND

Replaces RJ45 socket X6 of COMM-1

Screw terminals COMM2
COMM-H / COMM-L / GND

Replaces RJ45 socket X7 of COMM-2

Shielding contact area



Wiring notes:

- Uses either screw terminals or RJ45 sockets of one EASYPAB controller.
- Use always two cores for each signal
- Use clamp to contact cable shielding with shielding contact area on the controller PCB

Wire assignment as follows:

Screw terminal TCU3	Cable core (Color coding according EIA/TIA568B)
GND	blue and white/blue
COMM_H	white/orange and white/green
COMM_L	orange and green



Communication system - Communication line (CL) and it's termination

Important notes for system design



- Up to 24 EASYLAB controller can be connected together
 FH: Fume cupboard controller SC: Supply air controller EC: extract air controller RS: Room supply air controller
 RE: Room extract air controller, PC: pressure controller, TAM: TROX Adaption module

- Any order of devices at communication line is allowed

- Only controller types (RS, RE, PC, TAM) of **one** system type can work within one EASYLAB communication system; that means either controller types of an extract air controller system (LAB Laboratory) or controller types of a supply air controlled system (CLR CleanRoom)

- Maximum length of the communication line CL: 300m.

- Use network cable type SF-UTP acc. ISO IEC 11801 (braid+foil shielded e.g. CAT 5)

- 5 m connection cable for communication line CL included in the delivery
 (only for fume cupboard controller with device function FH-xxx)

- Communication line termination at the beginning and the end
 by activating the switch COMM-Terminator on the EASYLAB main board (see examples below)




Termination switch on the pcb

Symbol activated / deactivated


Example 1: Two fume cupboard controller or two room controller




Example 2: Two fume cupboard controller and one room supply air controller



Example 3: Three fume cupboard controller and one room supply air controller



Example 4: Two fume cupboard controller, one extract air controller, two room supply air controller and one room extract air controller





Room-Management-Function

The Room-Management-Function (RMF) is a software option to realize room based functions with the EASYLAB System:


- Centralized interchange point for wiring of room functions
- Connection of room control panels
- Centralized configuration of room parameters and functions
- Monitoring of room functions

Technical data:


- RMF can be activated on any room controller (Supply or Extract air) or on a TROX Adaption module.
- RMF activation is only allowed on exactly one EASYLAB controller in a system of max. 24 controller.
- RMF can be activated with the EASYLAB configuration software EasyConnect.



Example 1:




Example 2:




RMF relevance for wiring

The selection of the controller with RMF is relevant for the wiring of room based functions and signals. Following signals and the room control panel have to be connected at the controller with the RMF:





Connection of control panels



Note:

Use always blue patch cables for the control panels

Aim:

Simple determination of communication line cable (green) and control panel connecting cable (blue)



Control panel for fume cupboard controllers:



BE-SEG-02



CP-TOUCH

Technical data:

- Admitted are control panel types BE-SEG-01 and -02, BE-LCD and CP-TOUCH
- Up to two control panels can be connected to one controller
- 5 m connection cable blue included in the delivery
- Optionally connection cable up to 40 m overall length useable: Network cable type **SF-UTP** acc. ISO IEC 11801 (braid+foil shielded) with RJ45 connectors on both sides, e.g. patch cable Cat 5
- Plug in connection cable on back side of control panel

Room control panel:



CP-TOUCH

Technical data:

- Useable only at room controllers with activated Room-Management-Fuction
- Admitted is only control panel type BE-LCD and CP-TOUCH
- Up to two control panels can be connected to one controller
- 5 m connection cable blue included in the delivery
- Optionally connection cable up to 40 m overall length useable: Network cable type **SF-UTP** acc. ISO IEC 11801 (braid+foil shielded) with RJ45 connectors on both sides, e.g. patch cable Cat 5
- Plug in connection cable on back side of control panel

Note:

Room control panels can only be used at room controllers with activated Room-Management-Fuction (RMF).






Connection of TCU3 and PC for TROX EasyConnect-Software

EASYPANEL Controller / Control panel

PC-Software TROX EasyConnect

- Configuration
- Diagnostics
- Maintenance



EASYPANEL Configuration cable 5 m
TROX M-Number: M516SM3

USB-RS485 Adaption modul
TROX M-Number: M536ED7

BlueCON Adaption module
TROX M-Number: M546GA1

→ EASYPANEL Regler / Bedieneinheit

Complete package 1:
Order Code: **EasyConnect-CAB**

EasyConnect Software +
Configuration cable + USB-RS485 Adaption module +
programming adapter
TROX Number: B588NF4

Complete package 2:
Order Code: **EasyConnect-BC**

Alternatively you can use the Bluetooth Adaption module BlueCON for
a wireless connection between TCU3 and PC with EasyConnect
module and the cable.

EasyConnect Software + Bluetooth Adaption module BlueCON
+ programming adapter
TROX Number: B588NF5

Connection via control panels

BE-SEG-02



Programming adapter
TroX number: A00000043797


Control panel of fume cupboard controller

CP-TOUCH



Fume cupboard or room control panel

Connection directly at EASYPANEL controller



Use connection socket of control panels


RJ45 socket X3

RJ45 socket X2







Digital inputs (DI) / Digital outputs (DO)



Digital inputs DI1 ... DI6



Typical applications: Reading of switch contacts




The contact effect (Make/Brake contact) can be configured within the EasyConnect software


Connector for DI1:

- Phoenix MC 1,5/ 2-ST-3,5 Phoenix-No. 1840366 (included in the delivery)
TroX material number: A00000040301
- Phoenix MCVR 1,5/ 2-ST-3,5 Phoenix-No. 1863152 (by customer)
- Phoenix MCVW 1,5/ 2-ST-3,5 Phoenix-No. 1862852 (by customer)

Digital outputs DO1 ... DO6



Typical application: Special functions / Alarm relay



Technical data:


Electrical rating of contacts and pcb: 230 V AC max. 8 A

Additional note to alarm relay effect:


- Alarm relay released → Alarm state
- Alarm relay activated → No alarm state




Analogue inputs (AI) / Analogue outputs (AO) / Terminal 3




Analogue inputs AI1 ... AI5



Typical application: Sensor/Transducer integration



Typical application: Signal integration



Technical data:


- Input voltage range 0-10 V DC
- Characteristic freely configurable
- Input resistance > 100kOhm

Connector for AI5:


- Phoenix MC 1,5/ 3-ST-3,5 Phoenix-No. 1701879 (included in the delivery)
TroX material number: A00000040302
- Phoenix MCVR 1,5/ 3-ST-3,5 Phoenix-No. 1863165 (by customer)
- Phoenix MCVW 1,5/ 3-ST-3,5 Phoenix-No. 1862865 (by customer)



Analogue outputs AO1 ... AO4



Typical application: Signal transmission output



Technical data:

- Output voltage range 0-10 V DC
- Characteristic freely configurable
- Maximum output load 10 mA

Connector for AO4:

- Phoenix MVSTBR 2,5 /3-ST-5,08 Phoenix-No. 1792252 (included in the delivery)
- Trox material number: M516EE1

Terminal 3 (connection TROX HPD)




use: connection off digital actuator TROX HPD



Important note:

Actuator TROX HPD may be connected only to terminal 3!

Sash distance sensor DS-TRD





Input / Output assignment - Fume cupboard controller (FH)

Digital inputs DI		
DI 1	Sash monitoring EN14175	(only if configured; else unused)
DI 2	Switch contact of 2/3-point control	(only if device function is FH-2P or FH-3P, else unused)
DI 3	Switch contact of 3-point control	(only if device function is FH-3P, else unused)
DI 4	Fume scrubber operation request	(only if fume scrubber function configured; else unused)
DI 5	Support flow fan - Error	(only if supportive flow technology configured, else unused)
DI 6	not used	
<p>Unused inputs can be configured with the EasyConnect software for the usage of following functions:</p> <ul style="list-style-type: none"> - Integration of switchable constant volume flows into room balance - Integration of motion detector for special functions - Smoke extraction function - Operation mode preset High Mode, Low Mode, Shut off mode, Hand mode <p>(only available for fume cupboard controllers which do not follow the room operation mode due to special configuration) Operating mode changeover possible via control panel with configured keys The state of the not used inputs could be transmitted by the expansion module LON / BACnet / MODBUS (available from TCU3 Software Release 4.x)</p>		
Digital outputs DO		
DO 1	Alarm relay (1)	
DO 2	Fume scrubber – Enable (if configured, else unused)	
DO 3	Supportive flow technology fan - On/Off (if configured, else unused)	
DO 4	Automatic sash window closing - Open sash / Operation mode dependent switching e.g. fans (2) (if configured, else unused)	
DO 5	Automatic sash window closing - Close sash / Operation mode dependent switching e.g. fans (2) (if configured, else unused)	
DO 6	Fume cupboard interior light - On / Off (if configured, else unused)	
<p>(1) Alarm state = Relay released (2) From release 5.x : DO4 Standard, high, open, red. mode; b: DO4 Standard, high, open mode, DO5 red. mode Unused outputs relays could be energized via the expansion module LON / BACnet / MODBUS (available from TCU3 Software Release 5.x).</p>		
Analogue inputs AI		
AI 1	Internal volume flow transducer	(optionally external volume flow transducer)
AI 2	Temperature sensor (2)	(if temperature alarm is activated, else unused (1))
AI 3	unused (1)	
AI 4	Sash distance sensor	(only if device function is FH-DS or FH-DV, else unused)
AI 5	Face velocity transducer	(only if device function is FH-VS, else unused)
TI1	Temperature sensor PT1000 (2)	(if temperature alarm is activated)
<p>Unused inputs can be configured with EasyConnect software for integration of variable volume flow into room balance (1) from TCU V7.0 and EM-LON / EM-BAC-MOS-01 V4 or EM-IP V2 the voltage value could be transmitted via the communication card (2) from TCU V8.0</p>		
Analogue outputs AO		
AO 1	Current volume flow of controller	
AO 2	Total extract air room/Total supply air room/Set point supply air/voltage value via bus (1) (configurable)	
AO 3	Damper position of controller	
AO 4	Damper actuator	
(1) from Version 7.0		



Input / Output assignment – supply air / extract air control (SC/EC)

Digital inputs DI	
DI 1	not used, without function
DI 2	Switch contact of 2/3-point control <i>(only if device function is SC/EC-2P or SC/EC-3P, else without function)</i>
DI 3	Switch contact of 3-point control <i>(only if device function is SC/EC-3P, else without function)</i>
DI 4	Shut off mode <i>(configurable; else without function)</i>
DI 5	without function
DI 6	without function

The state of the inputs could be transmitted by the expansion module LON / BACnet / MODBUS (available from TCU3 Software Release 4.x).

Digital outputs DO	
DO 1	Alarm relay ⁽¹⁾
DO 2	Alarm relay, time dependent ⁽²⁾ <i>(configurable; else without function)</i>
DO 3	without function
DO 4	without function
DO 5	without function
DO 6	without function

⁽¹⁾ Alarm state = Relay released
⁽²⁾ Alarm relay time dependent: after release of an alarm the relay is energized for a settable time
Outputs without function could be energized via the expansion module LON / BACnet / MODBUS (available from TCU3 Software Release 5.x).

Analogue inputs AI	
AI 1	Internal volume flow transducer <i>(optionally external volume flow transducer)</i>
AI 2	without function ⁽¹⁾
AI 3	without function ⁽¹⁾
AI 4	setpoint volume flow <i>(only if device function is SC/EC E0/E2, else not used, without function)</i>
AI 5	not used, without function

⁽¹⁾ from TCU V7.0 and EM-LON / EM-BAC-MOS-01 V4 or EM-IP V2 the voltage value could be transmitted via the communication card

Analogue outputs AO	
AO 1	Current volume flow of controller
AO 2	Total extract air room/Total supply air room/Set point supply air/voltage value via bus ⁽¹⁾ <i>(configurable)</i>
AO 3	Damper position of controller
AO 4	Damper actuator

⁽¹⁾ from Version 7.0



Input / Output assignment - Room Controller (RS/RE/PC)

Note:

Some assignments are only available on room controllers with activated Room-Management-Function (RMF)

Digital inputs DI	
DI 1	unused
DI 2	unused <i>(Default at RMF: Room operation mode - Shut off mode)</i>
DI 3	unused <i>(Default at RMF: Room operation mode - Low mode)</i>
DI 4	unused <i>(Default at RMF: Room operation mode - High mode)</i>
DI 5	Pressure control - Door contact <i>(only at RMF, if internal pressure control is configured; else unused)</i>
DI 6	Pressure control - Changing setpoint <i>(only at RMF, if internal pressure control is configured; else unused)</i>
<p><i>Unused inputs can be configured with the EasyConnect software for the usage of following functions:</i></p> <ul style="list-style-type: none"> - <i>Integration of switchable constant volume flows into room balance</i> - <i>Operation mode preset for the room: High Mode, Low Mode, Shut off mode, Hand mode (only at RMF)</i> <p><i>The state of the not used inputs could be transmitted by the expansion module LON / BACnet / MODBUS (available from TCU3 Software Release 4.x).</i></p>	
Digital outputs DO	
DO 1	Alarm relay ⁽¹⁾
DO 2	Switchable output (On/Off) <i>(only at RMF; else unused)</i>
DO 3	Pressure alarm ⁽¹⁾ <i>(only at RMF, if internal pressure control is configured; else unused)</i>
DO 4	Summary alarm ⁽¹⁾ <i>(only at RMF; else unused)</i>
DO 5	Diffuser control / Sun blind – Open ⁽²⁾ / Room mode dependent ⁽³⁾ <i>(only at RMF, if configured; else unused)</i>
DO 6	Diffuser control / Sun blind – Close ⁽²⁾ / Room mode dependent ⁽³⁾ <i>(only at RMF, if configured; else unused)</i>
<p><i>⁽¹⁾ Relevant for alarm relay DO1, DO3, DO4: Alarm state = Relay released.</i></p> <p><i>⁽²⁾ From release 3.x; ⁽³⁾ From release 5.x, relay individual selectable</i></p> <p>Unused outputs relays could be energized via the expansion module LON / BACnet / MODBUS (available from TCU3 Software Release 5.x).</p>	
Analogue inputs AI	
AI 1	Internal volume flow transducer <i>(optionally external volume flow transducer)</i>
AI 2	unused ⁽¹⁾
AI 3	Volume flow offset - Temperature ⁽¹⁾ <i>(only at RMF, if function is configured; else unused)</i>
AI 4	Volume flow offset - Pressure <i>(only at RMF, if function is configured; else unused)</i>
AI 5	Press.Tranducer - Current pressure <i>(only at RMF, if internal pressure control is configured; else unused)</i>
<p><i>Unused inputs can be configured with EasyConnect software for integration of variable volume flow into room balance</i></p> <p><i>⁽¹⁾ from TCU V7.0 and EM-LON / EM-BAC-MOS-01 V4 or EM-IP V2 the voltage value could be transmitted via the communication card</i></p>	
Analogue outputs AO	
AO 1	Current volume flow of controller
AO 2	Total extract air room/Total supply air room/Set point supply air/voltage value via bus ⁽¹⁾ <i>(configurable)</i>
AO 3	Damper position of controller
AO 4	Damper actuator
<p>⁽¹⁾ from Version 7.0</p>	



Input / Output assignment - TROX Adaption Module (TAM)

Note:

Some assignments are only available on TAM's with activated Room-Management-Function (RMF)

Digital inputs DI	
DI 1	unused
DI 2	unused (Default at RMF: Room operation mode - Shut off mode)
DI 3	unused (Default at RMF: Room operation mode - Low mode)
DI 4	Unused (Default at RMF: Room operation mode - Hi mode)
DI 5	Pressure control - Door contact (only at RMF, if internal pressure control is configured; else unused)
DI 6	Pressure control - Changing setpoint (only at RMF, if internal pressure control is configured; else unused)

Unused inputs can be configured with the EasyConnect software for the usage of following functions:

- Integration of switchable constant volume flows into room balance
- Operation mode preset for the room: High Mode, Low Mode, Shut off mode, Hand mode (only at RMF)

The state of the not used inputs could be transmitted by the expansion module LON / BACnet / MODBUS (available from TCU3 Software Release 4.x).

Digital outputs DO	
DO 1	Alarm relay ⁽¹⁾
DO 2	Switchable output (On/Off) (only at RMF, else unused)
DO 3	Pressure alarm ⁽¹⁾ (only at RMF, if internal pressure control is configured, else unused)
DO 4	Summary alarm ⁽¹⁾ (only at RMF, else unused)
DO 5	Diffuser control / Sun blind – Open ⁽²⁾ / Room mode dependent ⁽³⁾ (only at RMF, if configured; else unused)
DO 6	Diffuser control / Sun blind – Close ⁽²⁾ / Room mode dependent ⁽³⁾ (only at RMF, if configured; else unused)

⁽¹⁾ Relevant for alarm relay DO1, DO3, DO4: Alarm state = Relay released.
⁽²⁾ From release 3.x; ⁽³⁾ From release 5.x, relay individual selectable
Unused outputs relays could be energized via the expansion module LON / BACnet / MODBUS (available from TCU3 Software Release 5.x).

Analogue inputs AI	
AI 1	unused
AI 2	unused ⁽¹⁾
AI 3	Volume flow offset - Temperature ⁽¹⁾ (only at RMF, if function is configured; else unused)
AI 4	Volume flow offset - Pressure (only at RMF, if function is configured; else unused)
AI 5	Press.Transducer - Current pressure (only at RMF, if internal pressure control is configured; else unused)

Unused inputs can be configured with EasyConnect software for integration of variable volume flow into room balance
⁽¹⁾ from TCU V7.0 and EM-LON / EM-BAC-MOS-01 V4 or EM-IP V2 the voltage value could be transmitted via the communication card

Analogue outputs AO	
AO 1	
AO 2	Total extract air room/Total supply air room/Set point supply air/voltage value via bus ⁽¹⁾ (configurable)
AO 3	
AO 4	

⁽¹⁾ from Version 7.0




Fume cupboard controller with face velocity transducer

Device function FH-VS

Technical data - transducer

- Supply voltage 24 V AC $\pm 15\%$ (via controller)
- Analogue output signal 2..10 V DC according to 0..1 m/s
- Connection cable approx. 3 m;
- Optionally, extension cable approx. 5 m; TROX Article-No. M536BA9

Face velocity transducer
Type VS-TRD



EN14175 Sash monitoring:
Bistable switch contact must be provided on site.

Fume cupboard controller with sash distance sensor


Device function FH-DS or FH-DV

Technical data - sensor

- Supply voltage 12 V DC (via controller)
- Analogue output signal according to extension length up to 1750 mm
- Connection cable approx. 3 m



Sash distance sensor
Type DS-TRD-02



EN14174 Sash monitoring
Monitoring of sash position
(500 mm) via sash distance
sensor.

Additional switch contact and
wiring not required!



Fume cupboard controller with face velocity transducer & sash distance sensor

Device function FH-VD


Technical data - transducer

- Supply voltage 24 V AC $\pm 15\%$ (via controller)
- Analogue output signal 2..10 V DC according to 0..1 m/s
- Connection cable approx. 3 m;
- Optionally, extension cable approx. 5 m; TROX Article-No. M536BA9



Technical data - sensor

- Supply voltage 12 V DC (via controller)
- Analogue output signal according to extension length up to 1750 mm
- Connection cable approx. 3 m



EN14174 Sash monitoring


Monitoring of sash position (500 mm) via sash distance sensor.

Additional switch contact and wiring not required!



Fume cupboard controller with 3 switching levels

Device function FH-3P




S1, S2 and EN14175 sash monitoring:
Bistable switch contacts must be provided on site.



Fume cupboard controller with 2 switching levels

Device function FH-2P




S1, S2 and EN14175 sash monitoring:
Bistable switch contacts must be provided on site.



Fume cupboard controller with fixed value control

Device function FH-F




EN14175 sash monitoring:
Bistable switch contact must be provided on site.



Special functions - Fume cupboard controller Part 1


Sash window monitoring EN 14175

Activation of function with EasyConnect configuration software




Activation of fume scrubber

Activation of function with EasyConnect configuration software




Control of supportive flow technology

Activation / configuration of function with EasyConnect configuration software



Control of automatic sash device


Activation of function with EasyConnect configuration software





Special functions - Fume cupboard controller Part 2

Operation mode dependent switching




2 variants possible:

a:
DO4 switches by Standard mode, High mode, Open mode, Low mode
DO5 not used

b:
DO4 switches by Standard mode, High mode, Open mode
DO5 switches by Low mode

Control of fume cupboard lighting



For safety reasons the TCU3 casing must be locked in such a way that it can only be opened with a tool when a 230 V power supply is connected.

A safety latch with a sticker is available, which can be ordered under EM-TRF-Zubehoer (A0000055540).

When delivered with transformer module EM-TRF or EM-TRF-USV, this is installed and is not required additionally.

Notes: Activation of function with EasyConnect configuration software

Control of fume cupboard lighting with expansion module EM-LIGHT

Expansion module EM-LIGHT ready to plug in of a 230 V AC fume cupboard lighting (no external wiring needed if transformer module is installed)




Mating connector not in scope of delivery of EM-LIGHT
type: WAGO WINSTA MIDI 770-113
Can be ordered under EM-CPL

Notes: Activation of function with EasyConnect configuration software. By delivered regulator with EM-LIGHT the function fume cupboard lighting is activated

Further details can be found in the installation and operating manual of the EM-LIGHT expansion module.

Activation of smoke extraction function




Notes: DIx: Function is available on all unused digital inputs DI1..DI6 by using Make or Brake switch contact.
Configuration of function with EasyConnect configuration software
Digital input, activation direction and EASYPYLAB Smoke extraction functionality (Damper open mode or damper shut off mode) must be configured




Special functions - Fume cupboard controller Part 3

Integration of constant volume flows into room balance




Notes: *DIx*: Function is available on all unused digital inputs DI1..DI6 by using Make or Brake switch contact.
Configuration of function settings with EasyConnect configuration software
Digital input, activation direction and required volume flow which should be integrated into the room balance during contact activation must be configured.

Integration of variable volume flows into room balance




Notes: *AIx*: Function is available on all unused analogue inputs AI by connecting an analogue signal with 0-10 V.
Configuration of function setting with EasyConnect configuration software
Analogue input, characteristic and required volume flow values (Extract or supply air) which should be integrated into the room balance must be configured.

Temperature alarm activ




Temperature alarm passiv



note: The function could be parametrized with the configuration software EasyConnect alternativ as an activ or a passiv signal.
The temperature alarm is available from Version 8.0 of EASYLAB

Activation of motion detector functions



Notes: *DIx*: Function is available on all unused digital inputs DI1..DI6 by using Make or Brake switch contact.
Configuration of function with EasyConnect configuration software
Digital input, activation direction and EASYLAB Motion detector functionality with parameters must be configured



Supply air- / Extract air controller 0-10 V / 2-10 V


Device function SC-E0; SC-E2; EC-E0; EC-E2

Setpoint volumeflow via 0 – 10 V (E0) or 2 – 10 V (E2)

Difference SC / EC:

SC: actual volume flow as supply via CL to room regulation system

EC: actual volume flow as exhaust via CL to room regulation system





Supply air- / Extract air controller 3-point-regulation -3P


Device function SC-3P; EC-3P

Setpoint volumeflow via 2 switch inputs, 3 setpoints

Difference SC / EC:

SC: actual volume flow as supply via CL to room regulation system

EC: actual volume flow as exhaust via CL to room regulation system





Supply air- / Extract air controller 2-point-regulation -2P


Device function SC-2P; EC-2P

Setpoint volumeflow via 1 switch inputs, 2 setpoints

Difference SC / EC:

SC: actual volume flow as supply via CL to room regulation system

EC: actual volume flow as exhaust via CL to room regulation system





Supply air- / Extract air controller Fixed volume flow -F


Device function SC-F; EC-F

Setpoint volumeflow fix parameterized

Difference SC / EC:

SC: actual volume flow as supply via CL to room regulation system


EC: actual volume flow as exhaust via CL to room regulation system






Special functions – Room controller Part 1

Control of on-site provided sun blind control device (only at room controller with activated RMF and alternative to Activation/Deactivation of supply air diffusers and Operation mode dependent switching)



Control of on-site provided lighting / device (only at room controller with activated RMF)




For safety reasons the TCU3 casing must be locked in such a way that it can only be opened with a tool when a 230 V power supply is connected.

A safety latch with a sticker is available, which can be ordered under EM-TRF-Zubehoer (A0000055540).


When delivered with transformer module EM-TRF or EM-TRF-USV, this is installed and is not required additionally.

Activation/Deactivation of supply air diffusers (only at supply air controller and alternative to sun blind control and Operation mode dependent switching)




Notes: Function is only available if sun blind control is not activated
Configuration of function settings with EasyConnect configuration software

Integration of constant volume flows into room balance



Notes: Dix: Function is available on all unused digital inputs DI1..DI6 by using Make or Brake switch contact.
Configuration of function settings with EasyConnect configuration software
Digital input, activation direction and required volume flow which should be integrated into the room balance during contact activation must be configured.

Integration of variable volume flows into room balance





Notes: Aix: Function is available on all unused analogue inputs AI by connecting an analogue signal with 0-10 V.
Configuration of function setting with EasyConnect configuration software
Analogue input, characteristic and required volume flow values (Extract or supply air) which should be integrated into the room balance must be configured.



Special functions – Room controller Part 2

Operation mode dependent switching

(alternative to sun blind control and Activation/Deactivation of supply air diffusers)



Individually adjustable for each relay
Relay switches in any given operating mode.
Combinations are possible




Internal pressure control

- Pressure control via supply air controller; (Extract air controlled system - Laboratory)
- Pressure control via extract air controller (Supply air controlled system - Clean room) available as of TCU3 controller with Software Version 3.0
- Current pressure signal integrateable via analogue signal 0-10V; Characteristic configurable
- 2 independant setpoits with related parameters configurable
- Setpoint changeover via switch contact at DI6 configurable alternatively via Expansion Modul via LonWorks®, BACnet MS/TP, Modbus-RTU, BACnet IP, Modbus IP
- Special functions via door switch contact at DI5 configurable (Alarmhandling, Ccontrol behaviouren)
- All connections and settings for pressure control within controller with Room-Management-Function (RMF)
- Seperate alarm output DO3 for pressure control loop


1. System design - Pressure control within extract air controlled room (e.g. Laboratory)

FH: Fume cupboard controller EC: extract air controller RE: Room extract air controller
PC: Executing pressure controller – Supply air

Example 1.1: RMF-Function activated on executing pressure controller




Example 1.2: RMF-Function activated on other room controller (not recommended)




2. System design - Pressure control within supply air controlled room (e.g. Clean room)

FH: Fume cupboard controller SC: Supply air controller RS: Room supply air controller
PC: Executing pressure controller – Extract air

Example 2.1: RMF-Function activated on executing pressure controller



Example 2.2: RMF-Function activated on other room controller (not recommended)




Recommendation:

In tight rooms and for optimal use of the diagnostic functions, the room management function (RMF) should be activated on the operating pressure controller (PC) .



Pressure control connections:



Power supply for pressure transducer

The power supply of the pressure transducer **can** be provided by the EASYLAB TCU3 via 24V connector of analogue input X5 AI Sensor.

Important:

- TCU3 supplied with 24 V AC → Connector 24V of X5 AI Sensor provides 24V AC
- TCU3 supplied with 24 V DC → Connector 24V of X5 AI Sensor provides 24V DC
- TCU3 supplied with 230V AC → Connector 24V of X5 AI Sensor provides 24V DC.

The above mentioned different supply voltage types must be observed depending on selected pressure transducer type.

Huba Pressure transducer type 699


This transducer is approved for AC and DC power supply type. Therefore no adaption is required.




Briem Pressure transducer type GB604 MF 0,75

This transducer is approved for AC and DC power supply type. Therefore no adaption is required.


Old Briem transducer must be manually adapted to the power supply type. Default delivery setting ist AC supply. Using a 230 V AC supply for the TCU3 via EM-TRF or EM-TRF-USV or 24 V DC supply implies a manual change of the supply to DC type. Therefore change jumper setting from AC to DC within the Briem transducer casing.



Tubing



LonWorks®-Interface (Expansion module EM-LON)




Technical data:

- LonWorks®-Interface FT10
- Network structure: Free Topology / Twisted Pair
- Simple network integration via internal connected double terminals
- Service-Key
- Status-LED for Service key, data transmission and data reception
- Mounting of expansion module within EASYPYLAB TCU3 controller box
- Supply voltage 5 V DC provided by EASYPYLAB controller TCU3
- Each network segment must be terminated with network terminator
- The maximum number of network devices is based on the LonWorks® specifications from Echelon. There are no additional restrictions by TROX.

Data interface: Standard Network Variable Types (SNVT)

- Different data points according to controller device type:
 - Room controller / TAM with activated Room-Management-Function (RMF) or Fume cupboard controller / Room controller / TAM without RMF
 - Usage on controller with activated RMF
 - > Access to datapoints of room
 - Usage on controller without RMF
 - > Access to datapoints of single controller

Connection assignment:



Expansion module built into EASYPYLAB TCU3 controller box

- 1 Expansion module EM-LON
- 2 LON-B
- 3 LON-A
- 4 Shielding
- 5 Network cables
- 6 LON-B
- 7 LON-A
- 8 Shielding

Network cable types:


Following cable types are recommended for LonWorks®-networks by Echelon:

- TIA 568A Category 5
- 8471 oder 85102 (Belden)
- Level IV cable
- JY(St)Y 2x2x0,8 (use only twisted wire pairs for LON-A and LON-B)

For more details see assembly and operating instruction manual of EM-LON.

The EM-LON extension module is only available as a spare part since 2022 (note availability of goods).

BACnet MS/TP interface (Expansion module EM-BAC-MOD-01) Modbus RTU interface



Technical data:

- BACnet MS/TP interface or Modbus RTU interface selectable via configuration switch
- Configurable network address and communication parameter
- Simple network integration via internal connected double terminals
- Status-LED for display of data transfer and communication error
- Mounting of expansion module within EASYPAB TCU3 controller box
- Supply voltage 5 V DC provided by EASYPAB controller TCU3
- Current regulations for network design and number of network device must be observed: BACnet Standard 135-2004 Modbus according EIA-485

In particular this includes:

- Network topology with simple line structure
- Usage of twisted pair copper cables with shielding
- Compliance of polarity A- and B+ for all network devices
- 120 Ohm resistor for network termination at first and last network device
- Network BIAS resistor for BACnet networks
- Max. 32 network devices in each network segment
- Configuration of individual network address for each device

Data interface:

- BACnet objects according PICS documentation or Modbus data registers
- Different data points according to controller device type: Room controller / TAM with activated Room-Management-Function (RMF) or Fume cupboard controller / Room controller / TAM without RMF


Usage on controller with activated RMF

-> Access to datapoints of room

Usage on controller without RMF

-> Access to datapoints of single controller

Connection assignment:



- 1 Expansion module EM-BAC-MOD-01
PCB illustration differs a little from original;
(Screw terminal representation is correct, see photo above)
- 2 B+ EIA-485
- 3 A- EIA-485
- 4 Shielding
- 5 Network cables
- 6 B+ EIA-485
- 7 A- EIA-485
- 8 Shielding

Expansion module built into EASYPAB TCU3 controller box

For more details see assembly and operating instruction manual of EM-BAC-MOD-01.



BACnet IP-interface (Expansion module EM-IP) Modbus IP-interface




Technical data:

- BACnet interface Schnittstelle IP or Modbus interface IP, switcheable via integral webserver
 - Configurable network address and transmission parameter
 - Two RJ45 10/100 Mbit Ethernet connections (up to five EM-IP modules can be linked by a daisy chain)
 - Status LED for indication of data transmission and errors
 - Mounting of expansion module within EASYLAB TCU3 controller box
 - Supply voltage 5 V DC provided by EASYLAB controller TCU3
- Data interface:

- BACnet objects according PICS documentation or Modbus data registers
 - Different data points according to controller device type:
Room controller / TAM with activated Room-Management-Function (RMF)
or Fume cupboard controller / Room controller / TAM without RMF
- Usage on controller with activated RMF
-> Access to datapoints of room
- Usage on controller without RMF
-> Access to datapoints of single controller

Connection assignment:



Expansion module EM-IP

- 1 Ethernet-socket LAN 1
- 2 Ethernet-socket LAN 2
- 3 LED-state indicator LAN 1
- 4 LED-state indicator LAN 2

Expansion module built into EASYLAB TCU3 controller box

For more details see assembly and operating instruction manual of EM-IP.